CGUARD™

Carotid Embolic Protective Stent System

Table of Contents

No.	Language	
1	English	2
2	Česky..	7
3	Deutsch.....	12
4	Dansk.....	18
5	Nederlands	23
6	Suomi.....	28
7	Français.....	33
8	Magyar	38
9	Italian.....	43
10	Latviski.....	48
11	Lietuviškai	53
12	Norsk.	58
13	Polish.....	63
14	Português.....	68
15	Русский язык.....	73
16	Srpski	79
17	Slovensky	84
18	Slovenščina	89
19	Español.....	94
20	Swedish	99


Figure 1 - CGuard Carotid Embolic Protection System (EPS)

Table of Contents

1	DEVICE DESCRIPTION	2
2	HOW SUPPLIED	2
3	INDICATIONS	2
4	CONTRAINDICATIONS	2
5	WARNINGS	3
5.1	General Warning	3
5.2	Device Use Warning	3
6	PRECAUTIONS	3
6.1	General –Precautions	3
6.2	Stent Handling – Precautions	3
6.3	Stent Placement – Precautions	3
6.4	Post-Implant – Precautions	4
7	MRI COMPATIBILITY	4
8	POTENTIAL ADVERSE EVENTS	4
9	STENT SIZE DETERMINATION	4
10	MATERIAL REQUIRED	4
11	PRE-PROCEDURE	5
12	INSPECTION PRIOR TO USE	5
13	DELIVERY SYSTEM PREPARATION	5
14	LESION PREPARATION	5
15	DELIVERY PROCEDURE	5
16	STENT DEPLOYMENT	5
17	POST STENT PLACEMENT	6
18	WARRANTY/LIABILITY	6

1 DEVICE DESCRIPTION

The InspireMD CGuard Carotid Embolic Protection System (EPS) is designed to deliver a self-expanding stent to the carotid arteries via an vascular sheath (7F, ID: 2.54 mm) or guiding catheter (8F, ID: 2.6mm), using an Over-the-Wire (OTW) Delivery System. The self-expanding stent is constructed of a nickel titanium alloy (Nitinol) and is covered by a permanent protective mesh. It is compressed and loaded into the OTW delivery system, delivered to the intended lesion site and then expanded by retraction of a protective sheath. The stent and mesh remain as a permanent vessel scaffolding implant. Upon deployment, the stent imparts an outward radial force on the arterial lumen to establish patency. The stents are available in the below size matrix.

Vessel diameter (mm)	Diameter (mm)	length (mm)			
		20	30	40	60
4.8 – 5.7	6.0	MCS0620	MCS0630	MCS0640	MCS0660
5.6 – 6.5	7.0	MCS0720	MCS0730	MCS0740	MCS0760
6.4 – 7.3	8.0	MCS0820	MCS0830	MCS0840	MCS0860
7.2 – 8.1	9.0	MCS0920	MCS0930	MCS0940	MCS0960
8.0 – 9.0	10.0	MCS1020	MCS1030	MCS1040	MCS1060

Table 1- CGuard Size Matrix

The stent delivery system is a coaxial Over-the-Wire system (OTW). The delivery system has a 120 cm working length and is compatible with 0.014" guidewires and embolic protection devices. Sheaths and guiding catheters with a minimum inner diameter of 2.54 mm are compatible. The delivery system is comprised of an inner shaft located inside the outer shaft lumen, which is the sliding tube component. Attached to the proximal end of the inner shaft is the luer, which permits flushing of the coaxial area and is the accesses of the guidewire. Two radiopaque marker bands are located on the inner shaft, one at each end of the stent to aid in positioning of the stent under fluoroscopy. A third radiopaque marker band located at the distal end of the outer sheath enables visualization of the sheath retraction during stent deployment. A pictorial representation of the Delivery System is presented in Figure 2.


Figure 2 - Delivery system components.
Above: Catheter and bifurcated luer,
middle: inner shaft and luer, below: inner tip

2 HOW SUPPLIED

This device is sterilized with ETO. Non-pyrogenic.
Contents: One (1) CGUARD Self-Expanding Carotid Stent with OTW Delivery System.
The system is placed inside a hoop and secured to the backing card for support.
Storage: Store in a dry, dark, cool place.

3 INDICATIONS

The CGuard Carotid Stent System is indicated for:
Improving carotid luminal diameter in patients at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet both criteria outlined below:

- Patients with neurological symptoms and >50% stenosis of the common or internal carotid artery by either ultrasound or angiogram OR patients without neurological symptoms and > 80% stenosis of the common or internal carotid artery by either ultrasound or angiogram, Patients having a vessel with reference diameters between 4.8 mm and 9.0 mm at the target lesion.

4 CONTRAINDICATIONS

The CGuard Carotid Stent System is contraindicated for use in:

- Patients in whom anti-coagulant and / or anti-platelet therapy is contraindicated.
- Patients with severe vascular tortuosity or anatomy that would preclude the safe introduction of a guide catheter, sheath, or stent system.
- Patients with known hypersensitivity to nickel-titanium.
- Patients with uncorrected bleeding disorders.
- Lesions in the ostium of the common carotid artery.
- Pre-existing jaw claudication.

5 WARNINGS

5.1 General Warning

- Only physicians who have received appropriate training and are familiar with the principles, clinical applications, complications, side effects and hazards commonly associated with carotid stent placement should use this device.
- As with any type of vascular implant, infection secondary to contamination of the stent may lead to thrombosis, pseudo-aneurysm, or rupture.
- Stenting across a major bifurcation may hinder or prevent future diagnostic or therapeutic procedures.
- When multiple stents are required, stent materials should be of similar composition.
- Do not use the product after the "Use By" date specified on the package.
- Do not use the product if the temperature indicator on inner pouch is black.
- Maintain the patient's Activated Clotting Time (ACT) at > 250 seconds throughout system usage to prevent thrombus formation on the device.
- Maintain continuous flush while removing and reinserting devices on the guidewire.
- Perform all exchanges slowly to prevent air embolism or trauma to the artery.
- Caution should be used if pre-dilating the lesion without embolic protection as this may increase the risk of an adverse outcome.
- Implanting a stent may lead to dissection of the vessel distal and / or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (carotid endarterectomy, further dilatation, or placement of additional stents).
- The stent may cause a thrombus, distal embolization or may migrate from the site of implant down the arterial lumen.
- Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infection, pseudo-aneurysm or fistulization, surgical removal of the stent may be required.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- If a filter (embolic protection system) is used, allow for and maintain adequate distance between the CGuard and the filter to avoid potential engagement or entanglement. If filter engagement and / or entanglement or filter detachment occurs, surgical conversion or additional catheter based intervention may be required.
- If a filter (embolic protection system) is used, choose a system with a long wire only (at least 300 cm).

5.2 Device Use Warning

- This device is intended for single-use only. DO NOT use the product after the 'Use By' date noted on the packaging.
- Package contains one self-expanding carotid stent system. Store at room temperature.
- DO NOT re-use. DO NOT re-sterilize, as this can compromise device performance and may increase the risk of cross-contamination due to inappropriate reprocessing.
- Ensure optimal positioning of the stent prior to deployment. Once deployment is, after initiated the stent cannot be repositioned or recaptured. Stent retrieval methods and the use of additional wires, snares, and/or forceps may result in additional trauma to the carotid vessel or the vascular access site. Complications may result in bleeding, hematoma, pseudoaneurysm, stroke or death.
- Overstretching of the artery may result in rupture and life-threatening bleeding.

- Implanting a stent may lead to dissection of the vessel distal to and/or proximal to the stent and may cause acute vessel closure requiring additional intervention such as carotid endarterectomy, placement of additional stents or further dilatation.
- Appropriate stent sizing is required to reduce the possibility of stent migration. The stent may migrate from the site of implant, embolize or cause a thrombus distally from the implant site down the arterial lumen.
- Maintain the patient's activated clotting time (ACT) at > 250 seconds throughout the stent delivery and implant procedure to prevent thrombus formation on the stent delivery system. In the event of thrombosis of the expanded stent, thrombolysis and PTA may be attempted. If IIb/IIIa inhibitors are used, maintain the ACT at > 200 seconds.
- Maintain a continuous flush while removing and inserting devices over the guidewire of the embolic protection device. Perform all device exchanges slowly in order to prevent air embolism or trauma to the artery.
- Allow for and maintain an adequate distance between the embolic protection system and the stent delivery system and/or the deployed stent in order to prevent possible entanglement of the two systems. If entanglement should occur and cannot be corrected, surgical intervention should be considered.

6. PRECAUTIONS

6.1 General- Precautions

Venous access should be available during carotid intervention to manage bradycardia or tachycardia by either medical therapy or temporary pacing.

6.2 Stent Handling - Precautions

- Special care must be taken not to handle or in any way disrupt the stent on the delivery system. This is most important during catheter removal from packaging, placement over the guidewire and advancement through the hemostasis valve and guiding catheter hub.
- Carefully inspect the System to verify that the device has not been damaged in shipment. Do not use damaged product.
- Take care to avoid unnecessary handling which may kink or damage the delivery system. Keep the delivery system as straight as possible and the delivery handle stationary during deployment. Do not use if device is kinked.
- Ensure the locking valve is open to assure freedom of movement of the outer sheath during deployment.
- Do not expose the delivery system to organic solvents as structural integrity and/ or device function may be impaired.
- Do not remove the stent from its delivery system, as removal may damage the stent. If removed, the stent cannot be put back onto the delivery system.
- The stent and the delivery system are designed to perform as an integrated system to be used only as designed.

6.3 Stent Placement - Precautions

- Ensure that the stent delivery system is fully flushed with heparinized saline prior to use. Do not use the delivery system if the flush is not observed exiting at the distal end of the catheter and through the locking valve.
- Use with bleedback control hemostatic valves is not recommended.
- The System is not compatible with guidewires or embolic protection devices larger than 0.014" wire (0.36 mm).
- The System must be used with a guiding catheter or sheath to maintain adequate support of the guidewire or embolic protection device throughout the procedure.

- If resistance is met during introduction or delivery of the stent delivery system, or during retraction of the outer sheath, the system should be carefully withdrawn and another system used.
- When catheters are in the body, they should be manipulated only under fluoroscopy. Radiographic equipment that provides high quality images is needed.
- Do not attempt to pull a partially deployed stent back through the guiding catheter or sheath; dislodgement of the stent from the delivery system may occur.
- Withdraw the stent system carefully if resistance is encountered during movement through the sheath or during initial retraction of the sheath.
- The delivery system is not designed for use with power injection. Use of power injection may adversely affect device performance.
- When more than one stent is required to cover the lesion or if there are multiple lesions, the distal lesion should be stented first followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent for placement of the distal stent and reduces the chance of dislodging stents that have already been placed.
- If overlap of sequential stents is necessary, the amount of overlap should be kept to a minimum (about 5 mm). In no instance should more than 2 stents overlap.

6.4 Post-Implant – Precautions

- Care must be exercised when crossing a newly deployed stent with other interventional devices to avoid disrupting the stent geometry and placement of the stent.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.

7 MRI COMPATIBILITY

The CGuard Self-expanding Nitinol stent system is MR conditional at field strength of 3 Tesla or less at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 min. Temperature rise produced should be less than 0.5°C.

8 POTENTIAL ADVERSE EVENTS

Based on the literature and on clinical and commercial experience with carotid stents and embolic protection systems, the following list includes possible adverse events associated with these devices:

- Abrupt closure
- Acute myocardial infarction
- Allergic reaction (contrast medium; drug; stent or filter material)
- Amaurosis fugax
- Aneurysm or pseudoaneurysm in vessel or at vascular access site
- Angina Coronary ischemia
- Arrhythmia (including premature beats, bradycardia, atrial and/or ventricular tachycardia, atrial and/or ventricular fibrillation [VF])
- Asystole or bradycardia requiring placement of a temporary pacemaker
- Arteriovenous fistula
- Bleeding complications from anticoagulant or antiplatelet medication requiring transfusion or surgical intervention
- Cerebral edema
- Cerebral hemorrhage
- Cerebral ischemia
- Congestive heart failure (CHF)
- Death
- Detachment and/or implantation of a component of the system
- Dissection of blood vessel
- Distal embolic protection device thrombosis occlusion

- Emboli, distal (air, tissue, plaque, thrombotic material, stent)
- Emergent or urgent surgery (Carotid Endarterectomy [CEA])
- Emergent surgery to remove stent or distal embolic protection device
- Fever
- Hematoma at vascular access site, with or without surgical repair
- Hemorrhagic event, with or without transfusion
- Hyperperfusion syndrome
- Hypotension/Hypertension
- Infection, local or systemic including bacteremia or septicemia
- Ischemia/ infarction of tissue organ
- Pain (head/neck)/ severe unilateral headache
- Pain at catheter insertion site
- Renal failure/insufficiency secondary to contrast medium
- Restenosis of vessel in stented segment
- Seizure
- Stent distal embolic protection device entanglement/ damage
- Stent distal embolic protection device collapse or fracture
- Stent malapposition/migration
- Stent thrombosis occlusion
- Stroke / cerebrovascular accident (CVA) / transient ischemic attack (TIA)
- Total occlusion of the carotid artery
- Vascular thrombosis/occlusion at puncture site, treatment site, or remote site
- Vessel dissection, perforation or rupture
- Vessel spasm or recoil

9 STENT SIZE DETERMINATION

Careful stent sizing is important to successful stenting. A minimum "interference" fit of 0.5 mm between the vessel and the stent is recommended in order to achieve optimum sizing and stent expansion of the self-expanding stent. For example, select a 6.0 mm stent to treat a 4.8 – 5.7 mm diameter vessel. Select a 7.0 mm stent to treat a 5.6 – 6.5 mm diameter vessel. The mean percentage of foreshortening for all stent sizes is less than 6%. The shortest stent length consistent with total lesion coverage is optimal. Should adequate coverage by one stent be impossible, a second stent may be used.

- The delivery system has a 120 cm working length and is compatible with 0.014" embolic protection devices and guidewires. Sheaths and guiding catheters with a minimum inner diameter of 2.54 mm are recommended.

WARNING: The CGUARD Carotid EPS is contraindicated for use with lesions in the ostium of the common carotid artery.

WARNING: Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration.

10 MATERIAL REQUIRED

- Vascular sheath (7F, ID: 2.54 mm) or guiding catheter (8F, ID: 2.6mm), guiding catheter or sheath length should not interfere with stent delivery system requirements.
- Optional balloon dilatation catheter.
- Any carotid EPD with a 300cm and 0.014" guidewire.
- Two to three syringes (10 -20 cc).
- 500 cc heparinized normal saline solution (sterile).

CAUTION: the CGUARD System is not compatible with guidewires or embolic protection devices larger than 0.014" (0.037 mm).

11 PRE-PROCEDURE

Patient preparation and sterile precautions should be the same as for any angioplasty procedure. The placement of the carotid stent in a stenotic or obstructed carotid artery must be done in a procedure room with angiography capabilities. Angiography should be performed to map out the extent of the lesion and the collateral flow. Access vessels must be sufficiently patent or sufficiently recanalized to proceed with further intervention.

12 INSPECTION PRIOR TO USE

- Remove the CGUARD Carotid EPS Self-Expanding Stent with OTW Delivery System from its protective packaging. Lay the device flat. Take care not to kink the shaft of the delivery catheter system.

CAUTION: Carefully inspect the CGUARD Carotid EPS Self-Expanding Stent with OTW delivery System to verify that the device has not been damaged in shipment. Do not use damaged equipment.

CAUTION: The delivery system has an internal shaft. Take care to avoid unnecessary handling which may kink or damage the delivery system. Keep the delivery system as straight as possible and the delivery handle stationary during deployment. Do not use if device is kinked.

- Inspect the delivery system sheath to verify that it has not been damaged during shipment and that the stent does not overlap the proximal marker.
- Ensure that the stent is fully covered by the sheath.

CAUTION: Special care must be taken not to handle or in any way disrupt the stent on the delivery system. This is especially important during delivery system removal from packaging, placement over the distal embolic protection device wire and advancement through a hemostatic valve and guiding catheter hub.

CAUTION: The stent on the delivery system is intended to perform as a system. Do not remove the stent from the delivery system as removal may damage the stent. If removed, the stent cannot be put back on the delivery system.

- Verify that the stent is the correct diameter and length by reading the specifications on the handle at the inner shaft luer of the delivery system. Do not use if any defects are noted.

13 DELIVERY SYSTEM PREPARATION

CAUTION: Do not expose the delivery system to organic solvents as structural integrity and / or function may be impaired.

- Fill a 10 cc syringe with heparinized normal saline and inject the saline through the luer fittings at the proximal end of the inner shaft luer for both the guidewire lumen and bifurcated luer at the device lumen. Flush until fluid is observed exiting the delivery system at the distal end of the catheter.

CAUTION: Ensure that the stent delivery system is fully flushed with heparinized saline prior to use. Do not use the delivery system if flush is not visible exiting at the distal end of the sheath.

- Keep the device straight and flat to avoid kinking the shaft.

14 LESION PREPARATION

WARNING: Administer heparin dose sufficient to maintain an ACT of > 250 secs to prevent thrombus formation on the devices.

CAUTION: Venous access should be available during carotid stenting to manage bradycardia and or hypotension by a pacemaker placement or pharmaceutical intervention, if needed.

CAUTION: The CGUARD EPS Self-Expanding Carotid Stent with OTW Delivery System must be used with a guiding catheter or sheath to maintain adequate support of the 0.014" guidewire

or embolic protection device throughout the procedure.

CAUTION: Use with bleedback control hemostatic valves is not recommended.

CAUTION: When catheters are in the, after body they should only be manipulated under fluoroscopy. Radiographic equipment that provides high-quality images is needed.

WARNING: Maintain continuous flush while removing and re-inserting devices on the guidewire or embolic protection device. Perform all exchanges slowly to prevent air embolism or trauma to the artery.

- If needed, pre-dilate the lesion with an appropriate size balloon dilatation catheter to a minimum of 3.0 mm after the distal protection device is in place beyond the lesion. *Note: If no predilatation balloon is utilized, there must be a minimum luminal opening of 3.0 mm to enable passage of the stent delivery system.*
- Maintain the embolic protection device wire and withdraw the balloon dilatation catheter.

15 DELIVERY PROCEDURE

- If lesion pre-dilatation has been performed, remove the balloon catheter and back-load the delivery system onto the 0.014" (0.36 mm) embolic protection device wire. A 300cm embolic protection device wire must be used. The wire will exit the delivery system at the luer.

CAUTION: The delivery system is not designed for use with a power injector. Use of a power injector may adversely affect device performance.

- Keep the device flat to avoid kinking the shaft.
- Insert the delivery system through the hemostatic valve adapter.

CAUTION: If resistance is met during delivery system introduction, the system should be withdrawn and another system used.

- Advance the stent and delivery system forward under fluoroscopic guidance to the lesion site.

CAUTION: Avoid tension in the delivery system prior to deployment.

16 STENT DEPLOYMENT

WARNING: Ensure optimal positioning of the stent prior to deployment. Once deployment is initiated, the stent cannot be repositioned or recaptured. Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the carotid vasculature and/or vascular access site. Complications may include death, stroke, bleeding, hematoma, or pseudoaneurysm.

CAUTION: Once stent placement has been initiated, do not attempt to pull a partially expanded stent back through the guiding catheter or sheath as dislodgement of the stent from the delivery system may occur.

- Confirm the stent position angiographically prior to deployment. Adjust position if necessary.
- While holding the handle stationary with one hand, unlock the locking valve.
- Note: Ensure that the delivery system is straight and not coiled. Keep the inner shaft of the delivery system stationary during deployment. Do not hold the outer sheath of the delivery catheter during deployment. It must be free to move.*
- Deploy the stent by slowly pulling back outer sheath. Continue to pull back until the stent is fully deployed.
- Note: If significant resistance is encountered during pull-back of the outer sheath and before stent release is initiated, re-lock valve and remove the system. Once deployment is initiated, the stent cannot be recovered by the sheath.*

CAUTION: In the event of partial delivery of the stent as the result of the inability to fully deploy the stent, remove the entire delivery system from the patient. This may result in damage to

the vessel wall and may require surgical intervention.

- Rotate the valve to the locked position, carefully withdraw the distal tip of the delivery system through the stent. Remove the delivery system from the patient.
- Under fluoroscopy, confirm that the stent has been deployed at the target lesion.
- If additional stent-to-wall apposition is desired or to facilitate the use of other interventional devices, the stent can be post-dilated with a balloon dilatation catheter. Do not expand the stent beyond its unconstrained maximum diameter as stated on the label and in Table 1. Post-dilate as needed in accordance with the compliance chart accompanying the selected balloon catheter.

CAUTION: When more than one stent is required to cover the lesion or if there are multiple lesions, the distal lesion should be stented first followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in order to place the distal stent and reduces the chance of dislodging stents that have already been placed.

CAUTION: If overlap of sequential stents is necessary, the amount of overlap should be kept to a minimum (approximately 5 mm). In no instance should more than two (2) stents overlap.

CAUTION: Care must be exercised when crossing a newly deployed stent with other interventional devices to avoid disrupting the stent geometry and placement of the stent.

WARNING: Overstretching of the artery may result in rupture and life-threatening bleeding.

17 POST STENT PLACEMENT

- Following stent placement, an angiogram should be performed to confirm vessel patency and percent stenosis remaining in the vessel lumen.


WARNING: The stent may cause a thrombus, distal embolization or may migrate from the site of implant down the arterial lumen. Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration. In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.

- Upon completion of the angiogram, the embolic protection device should be removed in accordance with the instructions for use with that device.
- Patients should be put on an appropriate regimen of anti-coagulants / antiplatelets.

WARNING: In the event of complications such as infection, pseudoaneurysm, or fistulization, surgical removal of the stent may be required.

18 WARRANTY/LIABILITY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in InspireMD instructions for use are expressly considered as an integral part of this provision. InspireMD warrants the product until the expiration date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. InspireMD disclaims any warranty of merchantability or fitness for a particular purpose of the product. InspireMD is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on InspireMD part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of InspireMD, however originating, from its product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects InspireMD legitimate interest in limiting its liability or warranty. No person has any authority to bind InspireMD to any warranty or liability regarding the product.

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InspireMD GmbH